

# Sports Medicine

CellRight Technologies®

As stewards of anatomical donations, CellRight Technologies® makes every effort to maximize this precious gift.



SPORTS MEDICINE GRAFTS	PRODUCT CODES
Whole Patellar Tendon	311001
Hemi-Patellar Tendon	311002
Pre-shaped Hemi-Patellar Tendon	311003
Achilles Tendon	311004
Anterior Tibialis Tendon	311006
Posterior Tibialis Tendon	311007
Semitendinosus	311008
Gracilis Tendon	311009
Peroneus	311010
Dermis	* call for sizes

Using our proprietary processing technology CellRight is able to increase the number of patients who can benefit from a transplant. Exceptional care and quality, coupled with scientific knowledge, ensures that we provide the foundation for the development of extraordinary regenerative products that positively impacts patient healing.

## MATRIX IQ® - DERMIS

Matrix IQ® human derived dermal grafts are decellularized to remove cellular components and preserve the biological properties that promote revascularization and repair of damaged tissue.

The indications for use of Matrix IQ® include the replacement of damaged or inadequate integument tissue, or for the repair, reinforcement or supplemental support of soft tissue defects such as breast reconstruction, abdominal wall repair, and reinforcement of tendons.

Matrix IQ® is approved as a wound covering and may be used in the treatment of diabetic ulcers, pressure ulcers, and partial and full-thickness wounds. Matrix IQ® maintains the elasticity of human skin and provides a natural scaffold for cellular and capillary remodeling of the damaged site.

## SAFETY

CellRight Technologies® employs strict quality assurance and quality

FRQWURO SURFHGXUHV WR HQVXUH SDWLHQW VDIHW\ &HOO5LJKW·V PHGL

director, a licensed physician, performs an extensive medical review of

WKH GRQRU·V PHGLFDO VRFLDO KLVWRU\ WR GHWHUPLQH HOLJLELWLW\ &H

Technologies® only accepts donors whose screening, serologic, and microbiologic tests meet or exceed the current standards established by the Food and Drug Administration (FDA) and the American Association of Tissue Banks (AATB).

Donor blood samples are tested by a CLIA licensed facility and must be found negative against infectious diseases such as Hepatitis B surface antigen, Hepatitis B core antibody, Hepatitis C antibody, Syphilis,

+,9 DQWLERG\ +,9 1\$7 +&9 1\$7 DQG +7/9 , , DQWLERG\

CellRight Technologies® products are processed aseptic in accordance with AATB standards and FDA federal regulations.

Grafts are sterilized using a validated low-dose, low-temperature gamma irradiation to achieve a medical device sterility assurance level (SAL) of 10<sup>-6</sup>. 6SRUWV PHGLFLQH JUDIWV DUH VWRUHG IURJHQ EHWZHHC

five-years from date of packaging.

